

LIPID MAPS Solution Protocol

Solution name: RAW 264.7 Serum Deprivation Medium 1

Solution name abbreviation: RAWSDM 1

Protocol ID: PS0000001000

Version: 00

Volume: 500ml

Components:

Reagent	Source	Catalog or Protocol No.	F.W. or Stock Conc.	Quantity	Final Concentration
DMEM with 4mM L-Glut and 4.5 g/L glucose	Cellgro	10-013	1X	497.5 ml	
Fetal Bovine Serum (FBS)	Hyclone	SH30071.03 ANG19242	100%	2.5 ml	.5%

To heat-inactivate the serum:

1. Thaw the 500 ml stock bottle of serum at 4°C.
2. Heat the serum at 56°C for 30 minutes.
3. Aliquot the serum in 50ml tubes and store at 4°C until use.

Solution Preparation:

1. Remove 2.5 ml from the 500ml DMEM reagent stock bottle.
2. Add 2.5 ml of heat-inactivated fetal calf serum.
3. Swirl to mix.
4. Warm the solution to 37°C before use.
5. Store at 4°C.

Author: Donna Reichart

Date: 11-7-05