

# Cellular Lipid Extraction for Targeted Stable Isotope Dilution Liquid Chromatography-Mass Spectrometry Analysis

Stacy L. Gelhaus, A. Clementina Mesaros, and Ian A. Blair. *J Vis Exp.* 2011; (57): 3399. Published online 2011 Nov 17. doi: 10.3791/3399

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3308612/>

## PROTOCOL

### 1. Standard and Internal Standard Mixes

1.1 Before lipid extraction, one must prepare a standard mix (SM) and an internal standard mix (ISM). For the SM, aliquot the equivalent volume of 1  $\mu\text{g}$  of each standard (25 total) listed in the reagent table below, using a calibrated syringe, into a 10 mL volumetric flask. Dry the standards under nitrogen and reconstitute in 10 mL of acetonitrile (ACN). The final concentration of the standard mix will be 100  $\text{pg}/\mu\text{L}$ . Divide into 1 mL aliquots and store at  $-80^{\circ}\text{C}$  until needed.

1.2 From this standard mix, 6 standard dilutions can be made for a total of 7 calibration standards (1  $\text{pg}/\mu\text{L}$ , 2  $\text{pg}/\mu\text{L}$ , 5  $\text{pg}/\mu\text{L}$ , 10  $\text{pg}/\mu\text{L}$ , 20  $\text{pg}/\mu\text{L}$ , 50  $\text{pg}/\mu\text{L}$ , and 100  $\text{pg}/\mu\text{L}$ ). These standards are to be run with each batch of samples prepared to develop the calibration curve used in quantitative analysis.

1.3 For the ISM, aliquot the equivalent volume of 1  $\mu\text{g}$  of each IS (14 total) listed in the reagent table below into a 10 mL volumetric flask, using a calibrated syringe, and dry down under nitrogen. Reconstitute the ISM in 10 mL of ACN and divide into 1 mL aliquots. The final concentration of the ISM should be 100  $\text{pg}/\mu\text{L}$ . There is enough ISM to process 1000 samples.

1.4 Store the SM and ISM at  $-80^{\circ}\text{C}$ .

### 2. Collection of Media and Cell Lysate

2.1 For each 10 cm plate of adherent cells ( $\sim 1 \times 10^7$ ), collect 3 mL of media in a 10 mL glass centrifuge tube. Next, wash the cells 2 times with phosphate buffered saline (PBS). Add 1 mL of PBS to the plate and scrape the cells off the plate and pipet into a second 10 mL glass centrifuge tube. At this point cells and media may be stored if the procedure will not be continued at this time. Store samples at  $-80^{\circ}\text{C}$  for long term storage or at  $-20^{\circ}\text{C}$  for less than one week.

### 3. Extraction of Lipids from Media

3.1 Free bioactive lipids are found in the cell media and these are easily extracted with a one-step l-l extraction procedure. To the collected media (2.1), add 10  $\mu\text{L}$  of ISM. At this time, 8 additional 10 mL glass centrifuge tubes should be labeled #0-7 for the calibration standards. Add 3 mL of PBS to each of these tubes and then add 10  $\mu\text{L}$  of the appropriate SM to tubes 1-7 (1-7 from 1.1).

3.2 Allow the ISM to equilibrate in each calibration standard and sample for 10 minutes and then add 5 mL of diethyl ether. The samples and calibration standards are then placed in the shaker on low speed for 30 minutes at room temperature.

3.3 After shaking, all calibration standards and samples are centrifuged at 1935 x g for 10 minutes. Next, using a Pasteur pipet, transfer the upper phase (organic) into a new 10 mL glass centrifuge tube and dry under nitrogen.

### 4. Extraction of Lipids from Cell Lysate

4.1 Spin down the cell samples in the glass centrifuge tubes at 1935 X g for 5 minutes. Dispose of the supernatant; add 2 mL of PBS to each tube, and vortex briefly. Remove 25  $\mu$ L of sample for protein quantification. Total protein will be used for cell lysate normalization.

4.2 Next, transfer half of the lysate/PBS mixture into a separate 10 mL glass centrifuge tube. Add an additional 2 mL of PBS along with 10  $\mu$ L of ISM to each sample. Follow the l-l extraction process described for media in 3.2 and 3.3. This aliquot will be used to measure free lipids in the cells and avoid alkaline degradation of lipids, such as the prostaglandins. This aliquot will not undergo steps 4.3 and 4.4 described below.

4.3 The first step of the esterified lipid l-l extraction process requires the addition of 5 mL of  $\text{CHCl}_3/\text{MeOH}$  (2:1) to each of the samples. Samples are set in a shaker on low speed for 30 minutes. After shaking, the samples are centrifuged at 1935 x g for 10 minutes. Next, using a Pasteur pipet, transfer the bottom (organic) layer of the l-l extraction into a new 10 mL glass centrifuge tube and dry the sample under nitrogen.

4.4 Esterified lipids must be saponified before quantification. Once the samples are dry, add 0.5 mL of 0.4N KOH in 80% MeOH to each tube. The samples are then incubated at 60°C for 1 h. After the incubation, add 2 mL of PBS. Adjust the pH to 6 with concentrated HCl acid (~30  $\mu$ L). Add 10  $\mu$ L of ISM to each sample. The rest of the procedure is the same as lipid extraction from media and can be followed in 3.2 and 3.3.

4.5 If these extractions are run simultaneously with the media samples the same calibration standards can be used to develop a calibration curve; however, if these samples are run at another time, calibration standards must also be processed as described in 3.1-3.3.

## 5. Pentafluorobenzylbromide (PFB) Derivatization

5.1 To the dried down samples and calibration standards, add 100  $\mu$ L of anhydrous dichloromethane ( $\text{CH}_2\text{Cl}_2$ ) followed by 100  $\mu$ L of 20% (v/v) N,N-diisopropylethylamine in anhydrous  $\text{CH}_2\text{Cl}_2$ , and finally 100  $\mu$ L of 20% (v/v) PFB in anhydrous  $\text{CH}_2\text{Cl}_2$ .

5.2 After the addition of derivatization reagents, the samples are vortexed, incubated at 60°C for 1 h and dried under nitrogen. The samples and calibration standards may be stored at -80°C if they are not analysed immediately or at -20°C for up to one week.

## 6. Stable Isotope Dilution Chiral Liquid Chromatography Electron Capture Atmospheric Pressure Chemical Ionization Mass Spectrometry (SID-LC-ECAPCI-MS) Analysis

6.1 The dried down samples and calibration standards are reconstituted in 100  $\mu$ L of hexane/ethanol (97:3). The entire sample can then be transferred to the HPLC vial with insert. Place the samples in the chilled autosampler (4°C) and create a sample list for analysis.

6.2 Next set up the HPLC and MS methods using the parameters described in Table 1 and Table 2 below. Table 3 contains specific compound information, including the MRM transition and collision energy (CE) for each lipid as well as the estimated retention time from the Waters 2695/Thermo Fisher TSQ Quantum Ultra analysis. The CE is determined by tuning each individual lipid during method development.

6.3 A Chiralpak AD-H column (4.6 X 250 mm, 5  $\mu$ ) heated to 30°C is used for the normal phase separation. The starting mobile phase is 98% hexane and 2% methanol/isopropanol (1:1). After the HPLC separation, methanol is added post-column (0.75 mL/min) to avoid deposition on the Corona needle that would cause a decrease in sensitivity. The MS is set in negative ion mode using the APCI source. If the source needs to be changed from ESI to APCI, be sure to insert a Corona needle.

6.4 The first sample that should be run is always a hexane blank. This will equilibrate the column and also reveal if there is any contamination on the column that may interfere with the calibration standard and sample analysis.

6.5 The next samples run should be the calibration standard samples in order of 0-7 (increasing SM concentration). These standards will be used to create a calibration curve containing the ratio of standard area/IS area on the y-axis and standard concentration or mass on the x-axis. From the  $y = mx + b$  equation derived from this plot, the ratio of sample area/IS area (y) can be used to determine sample (unknown) amount or concentration (x).



isoPs, and LTs. The collision energy (CE), MRM transitions, and retention time (RT, min) were optimized using a Waters 2695 LC coupled to a Thermo Fisher TSQ Quantum Ultra and are listed in Table 3 for each compound.

Table 1. HPLC gradient conditions for the chiral normal phase lipidomics separation. Solvent A is 100% hexanes and Solvent B is methanol/isopropanol (1:1).

<b>Time (min)</b>	<b>% A</b>	<b>% B</b>
0	98	2
3	98	2
11	96.4	3.6
15	92	8
27	92	8
30	50	50
38	50	50
39	98	2
45	98	2

Table 2. Listed above are the optimized mass spectrometry parameters for the measurement of bioactive lipids. These parameters are taken from a Thermo Fisher TSQ Quantum Ultra.

<b>Mass Spectrometry Parameter</b>	<b>Setting</b>
Source	APCI
Mode	Negative Ion
Method Length	40 min
MRM Scan Width	0.002 m/z
MRM Scan Time	0.15 sec
MRM Peak Width	0.70 FWHM
Tube Lens	164 V
Source CID	10 eV
Q2 Collision Gas	1.5 units
Vaporizer Temperature	450°C
Capillary Temperature	250°C
Discharge Current	0 $\mu$ A
Sheath Gas Pressure	25 units
Ion Sweep Pressure	3 units

Table 3. This table lists the compounds detected in the total ion chromatograms of Figure 1A-C. The second column provides the transition used for MRM/MS, the third column lists the collision energy (CE, eV) and the last column lists the retention time (RT, min). Parameters were optimized using a Waters 2695 separation module coupled to a Thermo Fisher TSQ Quantum Ultra.

Compound	Transition	CE (eV)	RT (min)
13-oxoODE	293.03 → 113.10	21	9.9
9-oxoODE	293.03 → 113.10	21	17.9
15-oxoETE	317.05 → 113.18	18	8.1
13(R)-HODE	295.03 → 195.12	18	12.5
13(S)-HODE	295.03 → 195.12	18	14.8
[2H4]-13(S)-HODE	299.03 → 198.12	18	14.9
9(R)-HODE	295.03 → 171.10	18	13.7
9(S)-HODE	295.03 → 171.10	18	21.2
[2H4]-9(S)-HODE	299.03 → 172.12	18	21.2
20-HETE	319.03 → 289.10	18	12.4
[2H6]-20-HETE	325.03 → 295.10	18	12.4
5(R)-HETE	319.04 → 115.04	16	15.7
5(S)-HETE	319.04 → 115.04	16	16.2
[2H8]-5(S)-HETE	327.04 → 116.04	16	16.4
12(R)-HETE	319.04 → 179.11	14	11.7
12(S)-HETE	319.04 → 179.11	14	12.3
[2H8]-12(S)-HETE	327.13 → 184.11	14	12.5
15(R)-HETE	319.04 → 219.12	13	12.4
15(S)-HETE	319.04 → 219.12	13	15.9
[2H8]-15(S)-HETE	327.03 → 226.12	13	16.1
11(R)-HETE	319.04 → 167.11	16	11.1
11(S)-HETE	319.04 → 167.11	16	12.5
8(R)-HETE	319.04 → 155.09	16	11.7
8(S)-HETE	319.04 → 155.09	16	12.6
LTB4	335.03 → 195.10	18	32
[2H4]-LTB4	339.03 → 197.10	18	32
PGE2	351.03 → 271.15	18	31.2
[2H4]-PGE2	355.03 → 275.15	18	31.2
PGD2	351.03 → 271.15	25	32.5
[2H4]-PGD2	355.03 → 275.15	25	32.5
8-iso-PGE2	351.03 → 271.15	18	37.4
11 $\beta$ -PGF2	353.03 → 193.11	25	27.6
[2H4]-11 $\beta$ -PGF2	357.03 → 197.11	25	27.6
PGF2 $\alpha$	353.03 → 193.11	25	33.1
[2H4]-PGF2 $\alpha$	357.03 → 197.11	25	33.1
8-iso-PGF2 $\alpha$	353.03 → 193.11	25	34.4
[2H4]-8-iso-PGF2 $\alpha$	357.03 → 197.11	25	34.4

## DISCUSSION

The standards and internal standards used in this protocol provide a representation of a targeted lipidomics method. A Waters 2695 separation module and Thermo Fisher TSQ Quantum Ultra were used for the LC-MS analysis and the optimal parameter settings can be found in Tables 1 and 2. Additionally, this extraction protocol was designed for adherent cells, but can be modified for other cell types as well as other biological matrices including urine, blood, and tissue. Many lipid standards are commercially available and can be incorporated into this analysis. For the derivatization procedure, it is important to use anhydrous solvents and the highest purity N,N-diisopropylethylamine. For the LC-MS analysis, it is also critical to use solvents of the highest purity. These solvents are listed in the reagent list. The post column addition of methanol is also critical to prevent deposition build up on the Corona needle as this would drastically decrease the sensitivity of the method.

The inclusion of calibration standards not only provides a means of accurate quantification, but also allows the researcher to troubleshoot problems with the extraction or LC-MS procedures. In some instances, a peak may be difficult to quantify. The sample chromatogram should be compared to the calibration standard chromatogram. It is often found that the calibration standard chromatogram does not contain this interference. This is indicative of an interfering compound, found in the cell media or lysate that is retained on the column at the same retention time as one of the lipids. The pH of the diethyl ether extraction in the protocol is 6. In some cases, optimizing the pH of this extraction step, such as extraction at acidic conditions (pH = 3) may remove interfering compounds. This parameter will be sample dependent as in many instances; the interfering species are specific to the cell type or medium used in the experiment.

Most current bioanalytical methods available for the measurement of bioactive lipids are not as extensive as they do not include chiral normal phase chromatography or SID. Chiral normal phase LC is critical for the separation of enantiomers and for being able to distinguish between enzyme- or ROS-mediated metabolism. The use of SID ensures that human error or complications that arise during extraction or analysis are taken into account during quantification. These added components along with ECAPCI-MRM make this the most sensitive, selective method available for the analysis of bioactive lipids.